

# OrthoBaby

MADE IN ITALY

A child's mouth is delicate  
Let's adapt to its nature



AA080444-ING

© Yareta



## New Asa Dental impression trays for children.

### **No plastic, all aluminium.**

A light, resistant, pliable and adaptable material.  
No more unpleasant burnt plastic smells or messy drilling problems.

**Low edges** to better adapt to the anatomy of little mouths.

**Reduced formats**, in 3 different sizes for school-age children.

Exactly what children want. And dentists, too.

**Children and dentists: smile!**

**Asa Dental S.p.A. uninominale**


Via Valenzana di Sopra, 60 - Bozzano - 55054 Massarosa (LU) Italy  
Tel. +39 0584 938305 - Fax +39 0584 937167 - [asadental@asadental.it](mailto:asadental@asadental.it)  
[www.asadental.it](http://www.asadental.it) - [www.asadental.com](http://www.asadental.com)

**ASA DENTAL**  
MAKE PEOPLE SMILE


# OrthoBaby

MADE IN ITALY


Aluminum impression trays  
for children



		A	B
		mm.	mm.
2816-U1	S	57	38
2816-U2	M	58	48
2816-U3	L	63	50



		A	B
		mm.	mm.
2816-L1	S	55	45
2816-L2	M	60	48
2816-L3	L	69	48



**S2816-6**  
MINI set of 6

Asa Dental S.p.A. uninominale

Via Valenzana di Sopra, 60 - Bozzano - 55054 Massarosa (LU) Italy  
Tel. +39 0584 938305 - Fax +39 0584 937167 - asadental@asadental.it  
www.asadental.it - www.asadental.com

**ASA DENTAL**  
MAKE PEOPLE SMILE